

Anti-Mouse Cd112/Nectin-2 Monoclonal Antibody (Abflo 488) For Fc - P32507

Item Number: CM0011382



Introduction

ABflo 488-conjugated rabbit anti-mouse CD112/Nectin-2 monoclonal antibody validated for flow cytometry. Ideal for T-cell signaling, immuno-oncology, and cell adhesion research in mouse models.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD112/Nectin-2 mAb
Remarks/Alias	MPH; Pvr; Pvs; Cd112; Pvr12; nectin-2
Species	Mouse
GeneID (human)	5819
GeneID	19294
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-345 of mouse CD112/Nectin-2 (NP_033016.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P32507
Protein Weight	57kDa
Shipping	Ice bag