

## Anti-Mouse Cd112/Nectin-2 Monoclonal Antibody For Fc - P32507

Item Number: CM0008364



### Introduction

Mouse CD112/Nectin-2 monoclonal antibody for flow cytometry, targeting the T-cell modulator and cell adhesion protein. Validated for mouse, ideal for immunology research.

[Learn More](#)

| Parameter        | Value                                                                                                                      |
|------------------|----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse CD112/Nectin-2 mAb                                                                                       |
| Remarks/Alias    | MPH, Pvr, Pvs, Cd112, Pvr12, nectin-2                                                                                      |
| Species          | Mouse                                                                                                                      |
| GeneID (human)   | 5819                                                                                                                       |
| GeneID (mouse)   | 19294                                                                                                                      |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 32-345 of mouse CD112/Nectin-2 (NP_033016.3) |
| Source           | Rabbit                                                                                                                     |
| Category         | Monoclonal Antibodies                                                                                                      |
| Application      | FC                                                                                                                         |
| Cross Reactivity | Mouse                                                                                                                      |
| SWISS            | P32507                                                                                                                     |
| Protein Weight   | 57kDa                                                                                                                      |
| Shipping         | Ice bag                                                                                                                    |