

Anti-Mouse Cd117/C-Kit Rabbit Monoclonal Antibody For Fc - P05532

Item Number: CM0009291



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD117/c-Kit receptor for flow cytometry. Supports hematopoietic and stem cell research, offered as a reliable IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD117/c-Kit mAb
Remarks/Alias	W; Bs; Fdc; Ssm; SC01; SC05; SOW3; CD117; c-KIT; Tr-kit; Gsfsc01; Gsfsc05; Gsfscow3
Species	Mouse
GeneID (Human)	3815
GeneID	16590
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 25-523 of mouse CD117/c-Kit (NP_001116205.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05532
Protein Weight	22kDa/108kDa/109kDa
Shipping	Ice bag