

Anti-Mouse Cd11A Monoclonal Antibody (Apc) For Fc - P24063

Item Number: CM0012810

Introduction

APC-conjugated rabbit monoclonal antibody targeting mouse CD11a (LFA-1 alpha chain) for flow cytometry. Recognizes leukocyte integrin subunit involved in adhesion and immune responses. Validated with recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD11a mAb
Remarks/Alias	Cd11a; LFA-1; Ly-15; Ly-21; (p180); LFA-1A
Species	Mouse
GeneID (Human)	3683
GeneID	16408
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 149-340 of mouse CD11a (NP_001240801.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P24063
Protein Weight	128kDa
Shipping	Ice bag