

Anti-Mouse Cd11A Monoclonal Antibody For Flow Cytometry - P24063

Item Number: CM0009265

Introduction

-conjugated rabbit anti-mouse CD11a monoclonal antibody for flow cytometry. Targets integrin alpha-L (LFA-1), involved in leukocyte adhesion, migration, and immune synapse formation. Suitable for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD11a mAb
Remarks/Alias	Cd11a; LFA-1; Ly-15; Ly-21; (p180); LFA-1A
Species	Mouse
Gene ID (Human)	3683
Gene ID (Mouse)	16408
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 149-340 of mouse CD11a (NP_001240801.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P24063
Protein Weight	128kDa
Shipping	Ice bag