


Anti-Mouse Cd11B Monoclonal Antibody (Biotin) For Flow Cytometry - P05555

Item Number: CM0004150



Introduction

Biotin-conjugated rabbit monoclonal antibody targeting mouse CD11b (Integrin alpha-M, ITGAM) for flow cytometry (FC) applications. Validated in mouse samples, recombinant immunogen, 127 kDa predicted weight.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Mouse CD11b mAb
Remarks/Alias	CR3, CR3A, MAC1, Cd11b, Ly-40, Mac-1, Mac-1a, CD11b/CD18, F730045J24Rik
Species	Mouse
GeneID (Human)	3684
GeneID	16409
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 143-337 of mouse CD11b (NP_032427.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05555
Protein Weight	127kDa
Shipping	Ice bag