


Anti-Mouse Cd11B Monoclonal Antibody () For Flow Cytometry - P05555

Item Number: CM0012723



Introduction

Rabbit anti-Mouse CD11b monoclonal antibody, conjugated to , for flow cytometry. Validated for mouse. Recombinant protein immunogen. Ideal for studying integrin alpha-M function in monocytes, macrophages, and microglia.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD11b mAb
Remarks/Alias	CR3; CR3A; MAC1; Cd11b; Ly-40; Mac-1; Mac-1a; CD11b/CD18; F730045J24Rik
Species	Mouse
GeneID (Human)	3684
GeneID (Mouse)	16409
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 143-337 of mouse CD11b (NP_032427.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05555
Protein Weight	127kDa
Shipping	Ice bag