


Anti-Mouse Cd11B Monoclonal Antibody For Flow Cytometry - P05555

Item Number: CM0005854



Introduction

Rabbit monoclonal antibody against mouse CD11b (Integrin alpha-M), validated for flow cytometry. Recognizes the extracellular domain. Suitable for studies of leukocyte adhesion, complement receptor function, and microglial activity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD11b mAb
Remarks/Alias	CR3; CR3A; MAC1; Cd11b; Ly-40; Mac-1; Mac-1a; CD11b/CD18; F730045J24Rik
Species	Mouse
Gene ID (Human)	3684
Gene ID	16409
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 143-337 of mouse CD11b (NP_032427.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05555
Protein Weight	127kDa
Shipping	Ice bag