


## Anti-Mouse Cd11B Monoclonal Antibody (Percp) For Flow Cytometry - P05555

Item Number: CM0006080



### Introduction

PerCP-conjugated rabbit monoclonal antibody against mouse CD11b, validated for flow cytometry. Targets integrin alpha-M (Mac-1), a leukocyte adhesion receptor. Suitable for immunological research. SWISS P05555.

[Learn More](#)

| Parameter        | Value                                                                   |
|------------------|-------------------------------------------------------------------------|
| Product Name     | PerCP Rabbit anti-Mouse CD11b mAb                                       |
| Remarks/Alias    | CR3; CR3A; MAC1; Cd11b; Ly-40; Mac-1; Mac-1a; CD11b/CD18; F730045J24Rik |
| Species          | Mouse                                                                   |
| GeneID (Human)   | 3684                                                                    |
| GeneID           | 16409                                                                   |
| Immunogen        | Recombinant protein                                                     |
| Source           | Rabbit                                                                  |
| Category         | Monoclonal Antibodies                                                   |
| Application      | FC                                                                      |
| Cross Reactivity | Mouse                                                                   |
| SWISS            | P05555                                                                  |
| Protein Weight   | 127kDa                                                                  |
| Shipping         | Ice bag                                                                 |