

Anti-Mouse Cd11C Mab For Fc - Q9Qxh4

Item Number: CM0012248

Introduction

Mouse CD11c monoclonal antibody conjugated to , validated for flow cytometry (FC). Suitable for studying integrin alpha-X in immune cell adhesion and inflammatory responses.

[Learn More](#)

Parameter	Details
Product Name	anti-Mouse CD11c mAb
Remarks/Alias	Cr4; N418; Cd11c
Species	Mouse
GeneID (Human)	3687
GeneID	16411
Immunogen	Recombinant protein of mouse CD11c
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QXH4
Protein Weight	129kDa
Shipping	Ice bag