

## Anti-Mouse Cd11C Monoclonal Antibody For Fc - Q9Qxh4

Item Number: CM0010837

### Introduction

anti-Mouse CD11c monoclonal antibody from Armenian Hamster, validated for Flow Cytometry (FC). Targets mouse Integrin alpha-X, a key adhesion receptor in immune responses. Ideal for immunology research.

[Learn More](#)

| Parameter        | Value                              |
|------------------|------------------------------------|
| Product Name     | anti-Mouse CD11c mAb               |
| Remarks/Alias    | Cr4; N418; Cd11c                   |
| Species          | Mouse                              |
| Gene ID (Human)  | 3687                               |
| Gene ID (Mouse)  | 16411                              |
| Immunogen        | Recombinant protein of mouse CD11c |
| Source           | Armenian Hamster                   |
| Category         | Monoclonal Antibodies              |
| Application      | FC                                 |
| Cross Reactivity | Mouse                              |
| SWISS            | Q9QXH4                             |
| Protein Weight   | 129 kDa                            |
| Shipping         | Ice bag                            |