

Anti-Mouse Cd11C Monoclonal Antibody For Fc - Q9Qxh4

Item Number: CM0008670



Introduction

Armenian hamster monoclonal antibody against mouse CD11c (Integrin alpha-X) validated for flow cytometry. Recognizes mouse CD11c; immunogen recombinant protein. Ideal for inflammatory response and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	anti-Mouse CD11c mAb
Remarks/Alias	Cr4; N418; Cd11c
Species	Mouse
GeneID (Human)	3687
GeneID	16411
Immunogen	Recombinant protein of mouse CD11c
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QXH4
Protein Weight	129kDa
Shipping	Ice bag