

Anti-Mouse Cd11C Monoclonal Antibody () For Flow Cytometry - Q9Qxh4

Item Number: CM0006015

Introduction

-conjugated anti-Mouse CD11c monoclonal antibody (clone N418) validated for flow cytometry. Targets integrin alpha-X, a dendritic cell marker, for mouse immune cell analysis.

[Learn More](#)

Parameter	Value
Product Name	anti-Mouse CD11c mAb
Remarks/Alias	Cr4; N418; Cd11c
Species	Mouse
GeneID (human)	3687
GeneID	16411
Immunogen	Recombinant protein of mouse CD11c
Source	Armenian Hamster
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QXH4
Protein Weight	129kDa
Shipping	Ice bag