

Anti-Mouse Cd124/Il-4Ra Monoclonal Antibody For Fc - P16382

Item Number: CM0009485

Introduction

Rabbit anti-Mouse CD124/IL-4Ra mAb validated for flow cytometry. Targets mouse interleukin-4 receptor alpha chain (UniProt P16382). Ideal for Th2 differentiation, allergy, and IgE studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD124/IL-4Ra mAb
Remarks/Alias	Il4r; CD124
Species	Mouse
Gene ID (Human)	3566
Gene ID	16190
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-233 of mouse CD124/IL-4Ra (NP_001008700.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P16382
Protein Weight	87kDa
Shipping	Ice bag