

Anti-Mouse Cd127/IL-7R α Monoclonal Antibody For Flow Cytometry - P16872

Item Number: CM0011468



Introduction

Rabbit monoclonal antibody targeting mouse CD127 (IL-7R α) for flow cytometry (FC). Validated in mouse, suitable for immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD127/IL-7R α mAb
Remarks/Alias	CD127; IL-7Ralpha
Species	Mouse
Gene ID (Human)	3575
Gene ID (Mouse)	16197
Immunogen	Recombinant protein (recombinant fusion protein containing a sequence corresponding to amino acids 21-239 of mouse CD127/IL-7R α , NP_612564.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P16872
Protein Weight	51kDa
Shipping	Ice bag