


Anti-Mouse Cd130/Gp130 Monoclonal Antibody For Fc - Q00560

Item Number: CM0008764



Introduction

-conjugated rabbit anti-mouse CD130/gp130 monoclonal antibody for flow cytometry (FC). Targets the IL-6 receptor subunit beta, enabling research on cytokine signaling, immune regulation, and inflammation. – Q00560.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD130/gp130 mAb
Remarks/Alias	CD130; gp130; D13Erttd699e; 5133400A03Rik
Species	Mouse
Gene ID (Human)	3572
Gene ID (Mouse)	16195
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-617 of mouse CD130/gp130 (NP_034690.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q00560
Protein Weight	102kDa
Shipping	Ice bag