


Anti-Mouse Cd130/Gp130 Monoclonal Antibody For Flow Cytometry - Q00560

Item Number: CM0008992



Introduction

Rabbit monoclonal antibody against mouse CD130/gp130 (SWISS Q00560), conjugated to dye, for flow cytometry. Reacts with mouse, immunogen maps to extracellular domain (aa 23-617).

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD130/gp130 mAb
Remarks/Alias	CD130; gp130; D13Ert699e; 5133400A03Rik
Species	Mouse
GeneID (Human)	3572
GeneID	16195
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-617 of mouse CD130/gp130 (NP_034690.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q00560
Protein Weight	102kDa
Shipping	Ice bag