

Anti-Mouse Cd138/Syndecan-1 Monoclonal Antibody () For Flow Cytometry - P18828

Item Number: CM0009776

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD138/Syndecan-1, validated for flow cytometry. Recognizes the extracellular domain, suitable for studying proteoglycan-mediated cell adhesion, exosome biology, and MSX1 signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD138/Syndecan-1 mAb
Remarks/Alias	Sstn; Synd; CD138; Synd1; syn-1
Species	Mouse
GeneID (Human)	6382
GeneID (Mouse)	20969
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-255 of mouse CD138/Syndecan-1 (NP_035649.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18828
Protein Weight	33kDa
Shipping	Ice bag