


Anti-Mouse Cd140B/Pdgfr Beta Monoclonal Antibody (Pe) For Flow Cytometry - P05622

Item Number: CM0008923



Introduction

PE-conjugated rabbit monoclonal antibody targeting mouse CD140b/PDGFR beta for flow cytometry. Suitable for angiogenesis, cardiovascular development, and renal fibrosis research. Validated in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Mouse CD140b/PDGFR beta mAb
Remarks/Alias	Pdgfr; CD140b; PDGFR-1
Species	Mouse
GeneID (Human)	5159
GeneID	18596
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-530 of mouse CD140b/PDGFR beta (NP_032835.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P05622
Protein Weight	123kDa
Shipping	Ice bag