

Anti-Mouse Cd146 Monoclonal Antibody () For Fc - Q8R2Y2

Item Number: CM0009184

Introduction

Rabbit anti-Mouse CD146 mAb validated for flow cytometry (FC). Targets mouse CD146/MUC18/gicerin, ~71kDa. Recombinant immunogen. Suitable for cell adhesion and melanoma research.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Mouse CD146 mAb
Remarks/Alias	CD146; CD149; Muc18; s-endo; 1-gicerin; s-gicerin
Species	Mouse
Gene ID (Human)	4162
Gene ID	84004
Immunogen	Recombinant protein corresponding to a sequence within amino acids 24-563 of mouse CD146 (NP_075548.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q8R2Y2
Protein Weight	71kDa
Shipping	Ice bag