

Anti-Mouse Cd150/Slam Monoclonal Antibody For Flow Cytometry - Q9Qum4

Item Number: CM0011136

Introduction

-conjugated rabbit anti-mouse CD150/SLAM monoclonal antibody for flow cytometry. Targets mouse SLAM family member 1, a 38 kDa cell membrane receptor involved in immune signaling. Reacts with mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD150/SLAM mAb
Remarks/Alias	Slam; CD150; IPO-3; CDw150; ESTM51; 4933415F16
Species	Mouse
GeneID (Human)	6504
GeneID	27218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-242 of mouse CD150/SLAM (NP_038758.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9QUM4
Protein Weight	38kDa
Shipping	Ice bag