

Anti-Mouse Cd19 Monoclonal Antibody () For Flow Cytometry - P25918

Item Number: CM0008152

Introduction

High-quality rabbit monoclonal antibody against mouse CD19, directly conjugated to for flow cytometry (FC). Validated in mouse samples; derived from recombinant mouse CD19 immunogen. Ideal for B-cell research applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD19 mAb
Species	Mouse
Gene ID (Human)	930
Gene ID	12478
Immunogen	Recombinant protein of mouse CD19.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25918
Protein Weight	60kDa
Shipping	Ice bag