


Anti-Mouse Cd204 Monoclonal Antibody For Flow Cytometry - P30204

Item Number: CM0005372



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD204 (MSR1), validated for flow cytometry. Ideal for macrophage scavenger receptor research in atherosclerosis, innate immunity, and foam cell formation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD204 mAb
Remarks/Alias	MSR; SR-A; Scvr; MRS-A; MSR-A; SR-AI; SR-AII; Scara1
Species	Mouse
Gene ID (Human)	4481
Gene ID	20288
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-458 of mouse CD204 (NP_001106797.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P30204
Protein Weight	50kDa
Shipping	Ice bag