

Anti-Mouse Cd215/Il-15R α Monoclonal Antibody () For If/Icc, Fc - Q60819

Item Number: CM0009236



Introduction

-conjugated rabbit monoclonal antibody specific for mouse interleukin-15 receptor subunit alpha (CD215/IL-15R α). Validated in IF/ICC and FC, it enables detection and analysis of IL-15 signaling in mouse samples. SWISS: Q60819, protein weight 28kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD215/IL-15R α mAb
Remarks/Alias	IL15ra
Species	Mouse
Gene ID (Human)	3601
Gene ID	16169
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 33-205 of mouse CD215/IL-15R α (NP_032384.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Mouse
SWISS	Q60819
Protein Weight	28kDa
Shipping	Ice bag