

Anti-Mouse Cd215/IL-15R α Monoclonal Antibody For If/Icc, Fc - Q60819

Item Number: CM0009560



Introduction

-conjugated rabbit anti-mouse CD215/IL-15R α monoclonal antibody for IF/ICC and flow cytometry. Specific to mouse; immunogen corresponds to aa33-205. Protein weight 28kDa. Suitable for IL-15 receptor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD215/IL-15R α mAb
Remarks/Alias	Il15ra
Species	Mouse
Gene ID (Human)	3601
Gene ID	16169
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Mouse
SWISS	Q60819
Protein Weight	28kDa
Shipping	Ice bag