

Anti-Mouse Cd22 Monoclonal Antibody () For Flow Cytometry - P35329

Item Number: CM0010842

Introduction

This -conjugated rabbit monoclonal antibody specifically detects mouse CD22 (Siglec-2), a B-cell inhibitory coreceptor, validated for flow cytometry in immunological and B-cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD22 mAb
Remarks/Alias	Lyb8; Lyb-8; A530093D23
Species	Mouse
GeneID (human)	933
GeneID	12483
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-151 of mouse CD22 (NP_033975.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35329
Protein Weight	89kDa/93kDa/97kDa
Shipping	Ice bag