

Anti-Mouse Cd22 Monoclonal Antibody () For If-P, Fc - P35329

Item Number: CM0008851

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD22 (Siglec-2). Validated for immunofluorescence and flow cytometry on B-cells, supporting B-cell biology and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD22 mAb
Remarks/Alias	Lyb8; Lyb-8; A530093D23
Species	Mouse
Gene ID (Human)	933
Gene ID (Mouse)	12483
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-151 of mouse CD22 (NP_033975.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, FC
Cross Reactivity	Mouse
SWISS	P35329
Protein Weight	89kDa/93kDa/97kDa
Shipping	Ice bag