

## Anti-Mouse Cd22 Rabbit Monoclonal Antibody For Flow Cytometry - P35329

Item Number: CM0011568

### Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD22 (Siglec-2). Validated for flow cytometry, ideal for B-cell biology, immunology, and signaling research. Detects native surface CD22.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse CD22 mAb                                                                                        |
| Remarks/Alias    | Lyb8, Lyb-8, A530093D23                                                                                           |
| Species          | Mouse                                                                                                             |
| GeneID (human)   | 933                                                                                                               |
| GeneID           | 12483                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 22-151 of mouse CD22 (NP_033975.3). |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | FC                                                                                                                |
| Cross Reactivity | Mouse                                                                                                             |
| SWISS            | P35329                                                                                                            |
| Protein Weight   | 89kDa/93kDa/97kDa                                                                                                 |
| Shipping         | Ice bag                                                                                                           |