

Anti-Mouse Cd223/Lag-3 Monoclonal Antibody For Fc - Q61790

Item Number: CM0012353



Introduction

-conjugated rabbit anti-mouse LAG-3 monoclonal antibody for flow cytometry. Targets mouse CD223, an immune checkpoint receptor on T-cells. Suitable for immunology research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD223/LAG-3 mAb
Remarks/Alias	Ly66; CD223; LAG-3
Species	Mouse
GeneID (Human)	3902
GeneID	16768
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-442 of mouse CD223/LAG-3 (NP_032505.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61790
Protein Weight	57 kDa
Shipping	Ice bag