

Anti-Mouse Cd23 Monoclonal Antibody For Fc - P20693

Item Number: CM0005061

Introduction

Rabbit anti-Mouse CD23 monoclonal antibody validated for flow cytometry (FC). Recognizes mouse CD23, low-affinity IgE receptor crucial for IgE regulation and B cell differentiation. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD23 mAb
Remarks/Alias	CD23; Fce2; Fcer2; Ly-42
Species	Mouse
GeneID (Human)	2208
GeneID	14128
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-331 of mouse CD23 (NP_038545.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P20693
Protein Weight	37kDa
Shipping	Ice bag