

Anti-Mouse Cd24 Monoclonal Antibody () For If/Icc, If-P, Fc - P24807

Item Number: CM0004092

Introduction

-conjugated rabbit anti-mouse CD24 monoclonal antibody validated for IF/ICC, IF-P, and flow cytometry. Suitable for mouse immune cell differentiation, thymocyte development, and B-cell activation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD24 mAb
Remarks/Alias	HSA; Cd24; Ly-52; nectadrin
Species	Mouse
Gene ID (Human)	100133941
Gene ID	12484
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-53 of mouse CD24 (NP_033976.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, IF-P, FC
Cross Reactivity	Mouse
SWISS	P24807
Protein Weight	8kDa
Shipping	Ice bag