

Anti-Mouse Cd24 Polyclonal Antibody For Wb, Fc, Elisa - P24807

Item Number: CM0023025

Introduction

mouse CD24 Rabbit pAb is a polyclonal antibody raised in rabbit against recombinant mouse CD24 protein. Validated for WB, FC, ELISA in mouse samples. Suitable for CD24 research in immunology, cell differentiation, and autoimmunity studies.

[Learn More](#)

Parameter	Value
Product Name	mouse CD24 Rabbit pAb
Species	Mouse
GeneID (Human)	100133941
GeneID	12484
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-53 of mouse CD24 NP_033976.1
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Mouse
SWISS	P24807
Protein Weight	8kDa
Shipping	Ice bag