

Anti-Mouse Cd252 (Ox40 Ligand) Monoclonal Antibody For If/Icc, Fc - P43488

Item Number: CM0005720

Introduction

High-quality rabbit monoclonal antibody targeting mouse CD252 (OX40 Ligand) for IF/ICC and flow cytometry. Ideal for T-cell co-stimulation research. Detects a 22kDa protein. UniProt P43488.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD252 (OX40 Ligand) mAb
Remarks/Alias	Ath1; gp34; Ath-1; Ox40l; TXGP1; CD134L; OX-40L; Tnlg2b; Txgp1l
Species	Mouse
GeneID (Human)	7292
GeneID	22164
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 51-198 of mouse CD252 (OX40 Ligand) (NP_033478.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Mouse
SWISS	P43488
Protein Weight	22kDa
Shipping	Ice bag