

Anti-Mouse Cd270/Hvem Monoclonal Antibody For Flow Cytometry - Q80Wm9

Item Number: CM0011009



Introduction

Rabbit anti-Mouse CD270/HVEM monoclonal antibody conjugated to , validated for Flow Cytometry. Suitable for studies of immune regulation, T cell signaling, and host defense. SWISS: Q80WM9.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Mouse CD270/HVEM mAb
Remarks/Alias	TR2; Atar; HveA; Hvem; Tnfrs14
Species	Mouse
GeneID (Human)	8764
GeneID	230979
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 39-207 of mouse CD270/HVEM (NP_849262.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q80WM9
Protein Weight	30kDa
Shipping	Ice bag