

Anti-Mouse Cd272/Btla Monoclonal Antibody () For Flow Cytometry - Q7Tsa3

Item Number: CM0010875



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD272/BTLA, validated for flow cytometry. Optimal for immunology and oncology research, detects BTLA on lymphocytes and other immune cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD272/BTLA mAb
Remarks/Alias	A630002H24
Species	Mouse
Gene ID (Human)	151888
Gene ID	208154
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-176 of mouse CD272/BTLA (NP_001032808.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q7TSA3
Protein Weight	34kDa
Shipping	Ice bag