


Anti-Mouse Cd274 (Pd-L1) Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Fc, Elisa - Q9Ep73

Item Number: CM0008156



Introduction

Rabbit monoclonal antibody against mouse CD274 (PD-L1, B7-H1). Validated for WB, IF/ICC, IHC-P, FC, ELISA. Targets immune checkpoint PD-L1 for cancer and autoimmune research. Mouse reactive.

[Learn More](#)

Parameter	Value
Product Name	CD274 Rabbit mAb
Remarks/Alias	B7h1; Pdl1; Pdcd1l1; Pdcd1lg1; A530045L16Rik
Species	Mouse
GeneID (human)	29126
GeneID	60533
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, FC, ELISA
Cross Reactivity	Mouse
SWISS	Q9EP73
Protein Weight	33kDa
Shipping	Ice bag