

Anti-Mouse Cd28 Monoclonal Antibody () For Flow Cytometry - P31041

Item Number: CM0012064



Introduction

Rabbit anti-Mouse CD28 monoclonal antibody conjugated to , validated for flow cytometry. Targets T-cell co-stimulatory receptor CD28 (Swiss-Prot P31041) for immunology research.

[Learn More](#)

Product Name	Rabbit anti-Mouse CD28 mAb
Remarks/Alias	Cd28
Species	Mouse
Gene ID (Human)	940
Gene ID (Mouse)	12487
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-149 of mouse CD28 (NP_031668.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P31041
Protein Weight	25kDa
Shipping	Ice bag