

Anti-Mouse Cd3 Monoclonal Antibody For Fc - P04235 / P22646

Item Number: CM0010092

Introduction

Rat anti-Mouse CD3 monoclonal antibody validated for flow cytometry, recognizing both CD3 delta and epsilon chains (P04235/P22646). Ideal for T-cell immunophenotyping and adaptive immunity research.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse CD3 mAb
Remarks/Alias	T3d; CD3; T3e; CD3epsilon
Species	Mouse
Gene ID (Human)	916/915
Gene ID	12500/12501
Immunogen	A synthesized peptide derived from mouse CD3
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P04235/P22646
Protein Weight	19kDa/21kDa
Shipping	Ice bag