

Anti-Mouse Cd3 Monoclonal Antibody For Flow Cytometry - P04235 / P22646

Item Number: CM0010174



Introduction

Rat anti-mouse CD3 monoclonal antibody validated for flow cytometry. Recognizes both CD3δ (19 kDa) and CD3ε (21 kDa) chains of the TCR-CD3 complex. Ideal for T-cell research in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rat anti-Mouse CD3 mAb
Remarks/Alias	T3d; CD3; T3e; CD3epsilon
Species	Mouse
GenelD	12500/12501
Immunogen	Synthetic peptide; A synthesized peptide derived from mouse CD3.
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P04235/P22646
Protein Weight	19 kDa/21 kDa
Shipping	Ice bag