


Anti-Mouse Cd317/Bst2 Monoclonal Antibody For Flow Cytometry - Q8R2Q8

Item Number: CM0008782



Introduction

-conjugated rabbit anti-mouse CD317/BST2 monoclonal antibody designed for flow cytometry applications. Specific for mouse CD317 (bone marrow stromal antigen 2), a key antiviral restriction factor expressed on type I interferon-producing cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD317/BST2 mAb
Remarks/Alias	GREG; Bst-2; CD317; DAMP-1; 2310015I10Rik
Species	Mouse
Gene ID (Human)	684
Gene ID	69550
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-152 of mouse CD317/BST2 (NP_932763.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Mouse
SWISS	Q8R2Q8
Protein Weight	19kDa
Shipping	Ice bag