


## Anti-Mouse Cd317/Bst2 Rabbit Monoclonal Antibody For Wb And Elisa - Q8R2Q8

Item Number: CM0005009



### Introduction

High-quality rabbit monoclonal antibody targeting mouse CD317/BST2 (bone marrow stromal antigen 2) for Western blot and ELISA. Suitable for antiviral and innate immunity studies in mouse models.

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD317/BST2 Rabbit mAb                                                                                                  |
| Remarks/Alias    | GREG; Bst-2; CD317; DAMP-1; 2310015I10Rik                                                                              |
| Species          | Mouse                                                                                                                  |
| GeneID (Human)   | 684                                                                                                                    |
| GeneID           | 69550                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 55-152 of mouse CD317/BST2 (NP_932763.1) |
| Source           | Rabbit                                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                                  |
| Application      | WB, ELISA                                                                                                              |
| Cross Reactivity | Mouse                                                                                                                  |
| SWISS            | Q8R2Q8                                                                                                                 |
| Protein Weight   | 19kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                |