

Anti-Mouse Cd32 Rabbit Monoclonal Antibody For Ihc-P, Fc, Elisa - P08101

Item Number: CM0012222

Introduction

High specificity anti-mouse CD32 rabbit monoclonal antibody validated for IHC-P, FC, ELISA. Suitable for studying Fc gamma receptor IIB functions in immune complex biology and B-cell modulation.

[Learn More](#)

Parameter	Value
Product Name	CD32 Rabbit mAb
Remarks/Alias	CD32; Fcgr2; Fcr-2; Fcr-3; Ly-17; LyM-1; Lym-1; fcRII; FcgRII; Fcgr2a; Ly-m20; Fc[g]RII; F630109E10Rik
Species	Mouse
GenelD (human)	2213
GenelD	14130
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-207 of mouse CD32 (NP_034317.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, FC, ELISA
Cross Reactivity	Mouse
SWISS	P08101
Protein Weight	31kDa/33kDa/36kDa
Shipping	Ice bag