

Anti-Mouse Cd34 Monoclonal Antibody For Flow Cytometry - Q64314

Item Number: CM0010426

Introduction

-conjugated rabbit monoclonal antibody against mouse CD34, a 41 kDa membrane protein. Validated for flow cytometry (FC), it detects hematopoietic progenitor cells. Ideal for hematopoiesis and stem cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD34 mAb
Remarks/Alias	CD34;CD34 molecule;GIG3;MORT1
Species	Mouse
GenelD (Human)	947
GenelD	12490
Immunogen	Recombinant Protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q64314
Protein Weight	41kDa
Shipping	Ice bag