


Anti-Mouse Cd38 Monoclonal Antibody For Fc - P56528

Item Number: CM0012440



Introduction

Rabbit monoclonal antibody against mouse CD38, validated for flow cytometry. Targets the ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 protein, a transmembrane glycoprotein involved in calcium signaling and lymphocyte activation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD38 monoclonal antibody
Remarks/Alias	I-19; ADPRC 1; Cd38-rs1
Species	Mouse
GeneID (Human)	952
GeneID	12494
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-304 of mouse CD38 (NP_031672.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P56528
Protein Weight	34kDa
Shipping	Ice bag