

## Anti-Mouse Cd4 Monoclonal Antibody (Apc/Cyanine7) For Flow Cytometry - P06332

Item Number: CM0009361

### Introduction

Rabbit monoclonal antibody against mouse CD4, conjugated to APC/Cyanine7, validated for flow cytometry. Ideal for studying T-helper cell development, immune response, and MHC class II coreceptor functions.

[Learn More](#)

| Parameter        | Value                                  |
|------------------|----------------------------------------|
| Product Name     | APC/Cyanine7 Rabbit anti-Mouse CD4 mAb |
| Remarks/Alias    | L3T4; Ly-4                             |
| Species          | Mouse                                  |
| GeneID (human)   | 920                                    |
| GeneID           | 12504                                  |
| Immunogen        | Recombinant protein                    |
| Source           | Rabbit                                 |
| Category         | Monoclonal Antibodies                  |
| Application      | FC                                     |
| Cross Reactivity | Mouse                                  |
| SWISS            | P06332                                 |
| Protein Weight   | 51kDa                                  |
| Shipping         | Ice bag                                |