

Anti-Mouse Cd4 Monoclonal Antibody (Pe/Cyanine5) For Fc - P06332

Item Number: CM0009903

Introduction

PE/Cyanine5 Rabbit anti-Mouse CD4 monoclonal antibody for flow cytometry (FC). Recognizes mouse CD4 (L3T4, Ly-4), a T-cell coreceptor for MHC class II. Validated for mouse samples, protein weight 51kDa. Suitable for immune response studies.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine5 Rabbit anti-Mouse CD4 mAb
Remarks/Alias	L3T4; Ly-4
Species	Mouse
GeneID (Human)	920
GeneID	12504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-394 of mouse CD4 (NP_038516.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06332
Protein Weight	51kDa
Shipping	Ice bag