

Anti-Mouse Cd45 Monoclonal Antibody () For Flow Cytometry - P06800

Item Number: CM0011631



Introduction

Rabbit anti-mouse CD45 monoclonal antibody (clone) validated for flow cytometry (FC). Detects leukocyte common antigen (L-CA, B220) in mouse samples. Ideal for immunophenotyping and T-cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD45 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
Gene ID (Human)	5788
Gene ID	19264
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 26-566 of mouse CD45 (NP_001104786.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130kDa / 140kDa / 144kDa
Shipping	Ice bag