

Anti-Mouse Cd45R/B220 Monoclonal Antibody (Biotin) For Flow Cytometry - P06800 / P08575

Item Number: CM0010933

Introduction

Biotin-conjugated rat monoclonal antibody against mouse CD45R/B220, validated for flow cytometry. Specific for mouse, recommended for immunophenotyping and T-cell signaling research. MW ~130-144 kDa.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rat anti-Mouse CD45R/B220 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyl-4
Species	Mouse
GeneID (Human)	5788
GeneID	19264/5788
Immunogen	Recombinant protein of mouse CD45R/B220
Source	Rat
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800/P08575
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag