

Anti-Mouse Cd5 Monoclonal Antibody (Apc) For Flow Cytometry - P13379

Item Number: CM0011488

Introduction

APC-conjugated rabbit monoclonal antibody against mouse CD5 (Ly-1/Lyt-1) for flow cytometry. Suitable for T-cell and B1a cell analysis, TCR/BCR signaling, and immune tolerance research.

[Learn More](#)

Parameter	Value
Product Name	APC Rabbit anti-Mouse CD5 mAb
Remarks/Alias	Ly-1; Ly-A; Ly-12; Lyt-1
Species	Mouse
GeneID (Human)	921
GeneID (Mouse)	12507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of mouse CD5 (NP_031676.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P13379
Protein Weight	54kDa
Shipping	Ice bag