

Anti-Mouse Cd50/Icam-3 (Icam-5, Telencephalin) Monoclonal Antibody () For Fc - Q60625

Item Number: CM0010435

Introduction

-conjugated rabbit monoclonal antibody against mouse CD50/ICAM-3 (ICAM-5, telencephalin) for flow cytometry. Targets neuronal adhesion molecule expressed in the brain; suitable for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD50/ICAM-3 mAb
Remarks/Alias	TLN; CD50; Tlcn; Icam3; ICAM-3; ICAM-5
Species	Mouse
GeneID (Human)	7087
GeneID (Mouse)	15898
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-227 of mouse CD50/ICAM-3 (NP_032345.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q60625
Protein Weight	97kDa
Shipping	Ice bag